

Decanoic Acid

Catalog No: #S3379

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Decanoic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	GluR antagonist
Calculated MW	?172.26
Formulation	C10H20O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:OthersSMILES:CCCCCCCCC(=O)OPathways:NeuroscienceReceptor:AMPABoling pt:Melting pt:Solubility:Soluble in Chloroform, DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pishan Chang, et al. Brain. 2016 Feb; 139(2): 431C443.

Note: This product is for in vitro research use only