

Lysionotin

Catalog No: #S3381


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lysionotin
Brief Description	Inhibitors
Purification	99.28%
Target Name	Others
Calculated MW	344.32
Formulation	C18H16O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:COC1=CC=C(C=C1)C2=CC(=O)C3=C(O2)C(=C(C(=C3O)OC)O)OCPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:Yellow powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Hu Y, et al. Mol Biol Rep. 2014 Mar;41(3):1693-702.

Note: This product is for in vitro research use only