

## Wilforlide A

Catalog No: #S3382


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Wilforlide A
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	454.7
Formulation	C30H46O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@]12CC[C@@H](C([C@H]1CC[C@@]3([C@H]2CC=C4[C@]3(CC[C@@]5([C@H]4C[C@]6[C[C@@H]5OC6=O)C)C)C)(C)C)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Xue M, et al. Fitoterapia. 2010 Dec;81(8):1109-12.

Note: This product is for in vitro research use only