

Cimifugin

Catalog No: #S3383


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cimifugin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	306.3
Formulation	C16H18O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC(C)([C@@H]1CC2=C(O1)C=C3C(=C2OC)C(=O)C=C(O3)CO)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etcAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Li Y, et al. Biomed Chromatogr. 2012 Oct;26(10):1234-40.

Note: This product is for in vitro research use only