


## Description

Product Name	TRX818
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	393.21
Formulation	C16H11FNNa2O5P
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:CancerSMILES:O=C1C=C(C2=CC=CC(F)=C2)NC3=C1C(P([O-])([O-])=O)=C(OC)C=C3.[Na+].[Na+]Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only