

Obacunone

Catalog No: #S3390


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Obacunone
Brief Description	Inhibitors
Purification	99.92%
Target Name	p38 MAPK inhibitor
Calculated MW	454.5
Formulation	C26H30O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@@]12CC[C@@H]3[C@]4(C=CC(=O)OC([C@@H]4CC(=O)[C@]3([C@@]15[C@H](O5)C(=O)O[C@H]2C6=COC(=C6)C)(C)C)CPathways:MAPKReceptor:p38 MAPKBoling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Murthy KN, et al.Toxicology. 2015 Mar 2;329:88-97.

Note: This product is for in vitro research use only