

Corosolic acid

Catalog No: #S3391


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Corosolic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	FAK inhibitor
Calculated MW	472.7
Formulation	C30H48O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@@H]1CC[C@@]2(CC[C@@]3(C(=CC[C@H]4[C@]3(CC[C@@H]5[C@@]4(C[C@H]([C@@H](C5(C)C)O)O)C)C)[C@@H]2[C@H]1C)C(=O)O

Pathways:Angiogenesis; Cytoskeletal SignalingReceptor:FAKBoling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Yoo KH, et al. Phytoter Res. 2015 May;29(5):714-23.

Note: This product is for in vitro research use only