

## Sinapine thiocyanate

Catalog No: #S3392


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sinapine thiocyanate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	368.45
Formulation	C16H24NO5CNS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C[N+](C)(C)CCOC(=O)/C=C/C1=CC(=C(C=C1OC)O)OC.C(#N)[S-] Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in Pyridine, Methanol, Ethanol, etc. Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Zhang QS, et al. Zhongguo Zhong Yao Za Zhi. 2014 Nov;39(22):4345-8.

Note: This product is for in vitro research use only