

## Tryptophenolide

Catalog No: #S3393


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tryptophenolide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	312.4
Formulation	C <sub>20</sub> H <sub>24</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others  
 SMILES: CC(C)C1=C(C2=C(C=C1)[C@]3(CCC4=C([C@@H]3CC2)COC4=O)C)O  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: 177-182°C  
 Solubility: Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.  
 Appearance: Cryst  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Rui Li, et al. International Journal of Mass Spectrometry, 2008, 278(1): 38C49.

Note: This product is for in vitro research use only