

Baohuoside I

Catalog No: #S3396


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Baohuoside I
Brief Description	Inhibitors
Purification	98.85%
Target Name	CXCR inhibitor
Calculated MW	514.52
Formulation	C27H30O10
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology
 SMILES: C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)OC2=C(OC3=C(C2=O)C(=CC(=C3CC=C(C)C)O)O)C4=CC=C(C=C4)OC)O)O
 Pathways: Autophagy; GPCR/G Protein; Immunology/Inflammation
 Receptor: CXCR4
 Boiling pt: Melting pt: Solubility: Soluble in Pyridine, Methanol, Ethanol, etc.
 Appearance: Yellow powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

- Kim B, et al. Biochemistry. 2014 Dec 9;53(48):7562-9.

Note: This product is for in vitro research use only