

Epimedin C

Catalog No: #S3397


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epimedin C
Brief Description	Inhibitors
Purification	99.11%
Target Name	Others
Calculated MW	822.8
Formulation	C39H50O19
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@@H]([C@H]([C@@H](O[C@H]2OC3=C(OC4=C(C3=O)C(=CC(=C4CC=C(C)C)O[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)O)C6=CC=C(C=C6)OC)C)O)O)OPathways:OthersReceptor:OthersBinding pt:Melting pt:Solubility:Pyridine, Methanol, Ethanol, etc.Appearance:Yellow powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Cui L, et al. Biomed Chromatogr. 2014 Oct;28(10):1306-12.

Note: This product is for in vitro research use only