

Ginsenoside Rg3

Catalog No: #S3402


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ginsenoside Rg3
Brief Description	Inhibitors
Purification	99.00%
Target Name	EGFR inhibitor
Calculated MW	785.02
Formulation	C42H72O13
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:CC(=CCC[C@@](C)([C@H]1CC[C@@]2[C@@H]1[C@@H](C[C@H]3[C@]2(CC[C@@H]4[C@@]3(CC[C@@H](C4(C)C)O[C@@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)[C@H]6[C@@H]([C@H]([C@@H]([C@H](O6)CO)O)O)C)O)C)O)CPathways:Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/AdaptorsReceptor:EGFRBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, Ethanol, etc.Appearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Joo EJ, et al. Chem Biol Interact. 2015 May 25;233:25-34.

Note: This product is for in vitro research use only