

6-Gingerol

Catalog No: #S3405


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	6-Gingerol
Brief Description	Inhibitors
Purification	99.31%
Target Name	Others
Calculated MW	294.4
Formulation	C17H26O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CCCC[C@@H](CC(=O)CCC1=CC(=C(C=C1)O)OC)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tzeng TF, et al. Phytomedicine. 2015 Apr 15;22(4):452-61.

Note: This product is for in vitro research use only