

Rapastinel

Catalog No: #S3407


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rapastinel
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others agonist
Calculated MW	413.47
Formulation	C18H31N5O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:C[C@@H](O)[C@@H](C(N)=O)NC([C@H]1N(C([C@H]2N(C([C@H]([C@H](O)C)N)=O)CCC2)=O)CCC1)=O
 Pathways:OthersRecept
 or:OthersBolng pt:Melting pt:Solubility:DMSO: 32 mg/MIAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and
 shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Burgdorf J, et al. Neuroscience. 2015 Nov 12;308:202-11.

Note: This product is for in vitro research use only