

Syringin

Catalog No: #S3408


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Syringin
Brief Description	Inhibitors
Purification	99.00%
Target Name	NOS inhibitor
Calculated MW	372.4
Formulation	C17H24O9
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: COC1=CC(=CC(=C1O[C@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)OC)/C=C/COP

Pathways: Immunology/Inflammation Receptor: NOS

Boling pt: Melting pt: Solubility: Soluble in Pyridine, Methanol, Ethanol, etc.

Appearance: Powder

Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Cui Y, et al. Fundam Clin Pharmacol. 2015 Apr;29(2):178-84.

Note: This product is for in vitro research use only