

Momordin Ic

Catalog No: #S3410


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Momordin Ic
Brief Description	Inhibitors
Purification	98.00%
Target Name	MAPK; PI3K
Calculated MW	764.93
Formulation	C41H64O13
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@]12CC[C@@H](C([C@H]1CC[C@@]3([C@H]2CC=C4[C@]3(CC[C@@]5([C@H]4CC(CC5)(C)C)C(=O)O)C)C)O[C@H]6[C@@H]([C@H]([C@H]([C@H](O6)C(=O)O)O)O[C@H]7[C@@H]([C@H]([C@H]([C@@H](CO7)O)O)O)OPathways:MAPK; PI3K/Akt/mTOR signalingReceptor:MAPK; PI3KBoling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, Ethanol, etc.Appearance:Brown powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang J, et al. Apoptosis. 2013 Jun;18(6):751-65.

Note: This product is for in vitro research use only