

Curdione

Catalog No: #S3411


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Curdione
Brief Description	Inhibitors
Purification	99.00%
Target Name	P450 inhibitor
Calculated MW	236.34
Formulation	C15H24O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1=CC(=O)C(C)C2=CC(=O)C(C)C12 Pathways: Metabolism Receptor: CYP3A4 Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Li J, et al. Asian Pac J Y Prev. 2014;15(22):9997-10001.

Note: This product is for in vitro research use only