

3,4-Dihydroxyphenylethanol


Catalog No: #S3412

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3,4-Dihydroxyphenylethanol
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	154.16
Formulation	C ₈ H ₁₀ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC(=C(C=C1CCO)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Brown Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Yunden J, et al. Chem Res Toxicol. 2009 May; 22(5): 835C841.

Note: This product is for in vitro research use only