

Morroniside

Catalog No: #S3414


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Morroniside
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others inhibitor
Calculated MW	406.38
Formulation	C17H26O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1[C@@H]2[C@H](CC(O1)O)C(=CO[C@H]2O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)C(=O)OCPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Xiong S, et al. J Pharmacol Sci. 2015 Apr;127(4):414-8.

Note: This product is for in vitro research use only