

Secoisolariciresinol diglucoside

Catalog No: #S3415


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Secoisolariciresinol diglucoside
Brief Description	Inhibitors
Purification	99.00%
Target Name	ROS inhibitor
Calculated MW	686.71
Formulation	C32H46O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

systemSMILES:COC1=C(C=CC(=C1)C[C@@H](CO[C@@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)[C@@H](CC3=CC(=C(C=C3)O)OC)O[C@@H]4[C@@H]([C@H]([C@@H]([C@H]([C@H](O4)CO)O)O)O)OPathways:Immunology/InflammationReceptor:ROSBoling pt:Melting pt:Solubility:Pyridine, Methanol, EthanolAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tao YL, et al. Appl Microbiol Biotechnol. 2014 Mar;98(6):2519-31.

Note: This product is for in vitro research use only