

Aucubin

Catalog No: #S3416


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aucubin
Brief Description	Inhibitors
Purification	99.00%
Target Name	IL Receptor; NF-κB
Calculated MW	346.33
Formulation	C15H22O9
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CO[C@H]([C@H]2[C@@H]1[C@@H](C=C2CO)O)[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O

Pathways: Immunology/Inflammation; NF-κB Receptor: IL-1β; NF-κB

Boling pt: Melting pt: Solubility: DMSO: 31 mg/mL Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wang SN, et al. Int Immunopharmacol. 2015 Feb;24(2):408-15.

Note: This product is for in vitro research use only