

Integrin α 3 Polyclonal Antibody

Catalog No: #41067

Package Size: #41067-1 50ul #41067-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin α 3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Integrin α 3 Polyclonal Antibody detects endogenous levels of Integrin α 3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Integrin α 3.
Target Name	Integrin α 3
Other Names	ITGA3; MSK18; Integrin alpha-3; CD49 antigen-like family member C; FRP-2; Galactoprotein B3; GAPB3; VLA-3 subunit alpha; CD antigen CD49c
Accession No.	Swiss-Prot: P26006NCBI Gene ID: 3675
Uniprot	P26006
GeneID	3675;
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of HepG2 cells using Integrin α 3 Polyclonal Antibody



Western Blot analysis of 293 HepG2 RAT-STOMACH cells using Integrin $\alpha 3$ Polyclonal Antibody

Note: This product is for in vitro research use only