

Ezutromid

Catalog No: #S3424


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ezutromid
Brief Description	Inhibitors
Purification	98.47%
Target Name	Others
Calculated MW	337.39
Formulation	C19H15NO3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:O=S(C1=CC=C(OC(C2=CC=C3C=CC=C3=C2)=N4)C4=C1)(CC)=O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mg/MIAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Tinsley JM et al. PLoS One. 2011 May 6;6(5):e19189.

Note: This product is for in vitro research use only