

## Polyphyllin?VII

Catalog No: #S3425


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Polyphyllin?VII
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	1049.21
Formulation	C51H84O22
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@H]1[C@H]2[C@H](C[C@@H]3[C@@]2(CC[C@H]4[C@H]3CC=C5[C@@]4(CC[C@@H](C5)O[C@@H]6[C@H]([C@@H]([C@H]([C@@H](O6)CO)O[C@H]7[C@@H]([C@H]([C@H](O7)CO)O)O[C@H]8[C@@H]([C@@H]([C@H]([C@H](O8)C)O)O)O)C)O[C@]1(CC[C@@H](C)CO[C@H]9[C@H]([C@H]([C@@H]([C@H](O9)CO)O)O)O)OCPPathways:OthersReceptor:OthersBoling pt:Melting pt:Solubility:Soluble in MethanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Li ZH, et al. Parasitology. 2013 Jul;140(8):952-8.;2. Yin X, et al.Biomed Chromatogr. 2013 Mar;27(3):343-8.;

Note: This product is for in vitro research use only