

Polydatin

Catalog No: #S3427


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Polydatin
Brief Description	Inhibitors
Purification	99.00%
Target Name	NF-κB inhibitor
Calculated MW	390.39
Formulation	C ₂₀ H ₂₂ O ₈
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: C1=CC(=CC=C1/C=C/C2=CC(=CC(=C2)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)O Pathways: NF-κB Receptor: NF-κB Boiling pt: Melting pt: Solubility: DMSO: 29 mg/mL Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ravagnan G, et al. Inflammation. 2013 Feb;36(1):26-34.

Note: This product is for in vitro research use only