

## Helicid

Catalog No: #S3428


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Helicid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	284.26
Formulation	C13H16O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C1=CC(=CC=C1C=O)O[C@H]2[C@@H]([C@@H]([C@@H]([C@H](O2)CO)O)O)OPathways:OthersReceptor:PXRBoling  
 pt:Melting pt:Solubility:Soluble in DMSOAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Chen Q, et al. Evid Based Complement Alternat Med. 2015, (2015): 9.

Note: This product is for in vitro research use only