

Cinobufagin

Catalog No: #S3432


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cinobufagin
Brief Description	Inhibitors
Purification	99.00%
Target Name	ATPase inhibitor
Calculated MW	442.55
Formulation	C26H34O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cancer/

ImmunologySMILES:CC(=O)O[C@@H]1[C@@H]([C@]2(CC[C@H]3[C@H]([C@@]24[C@@H]1O4)CC[C@H]5[C@@]3(CC[C@@H](C5)O)C)C6=

COC(=O)C=C6Pathways:Membrane transporter/Ion channelReceptor:Na+/K+-ATPaseBoiling pt:Melting pt:Solubility:Soluble in

DMSOAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Toma S, et al. Xenobiotica. 1987 Oct;17(10):1195-202.

Note: This product is for in vitro research use only