


Description

Product Name	TUG-891
Brief Description	Inhibitors
Purification	99.21%
Target Name	GPR agonist
Calculated MW	364.41
Formulation	C23H21FO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC1=CC=C(C2=CC=C(F)C=C2COC3=CC=C(CCC(O)=O)C=C3)C=C1 Pathways: Endocrinology/Hormones Receptor: GPR120 Bolin
 g pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Hopkins MM, Zhang Z, Liu Z, Meier KE. J Clin Med. 2016 Jan 26;5(2).

Note: This product is for in vitro research use only