

LYS01 free base

Catalog No: #S3436


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	LYS01 free base
Brief Description	Inhibitors
Purification	99.35%
Target Name	Others inhibitor
Calculated MW	440.368
Formulation	C23H23Cl2N5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CN(CCNC1=CC=NC2=CC(Cl)=CC=C12)CCNC3=CC=NC4=CC(Cl)=CC=C34
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Amaravadi RK, Winkler JD. Lys05: a new lysosomal autophagy inhibitor. *Autophagy*. 2012 Sep;8(9):1383-4.

Note: This product is for in vitro research use only