


Description

Product Name	S107
Brief Description	Inhibitors
Purification	99.63%
Target Name	Others
Calculated MW	209.31
Formulation	C11H15NOS
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CN1CCSC2=C(C1)C=C(C=C2)OCPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 38 mg/MIAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lehnart SE, et al. J Clin Invest. 2008 Jun;118(6):2230-45.

Note: This product is for in vitro research use only