

2-PMPA

Catalog No: #S3440


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-PMPA
Brief Description	Inhibitors
Purification	98.00%
Target Name	Carboxypeptidase(CPG2) inhibitor
Calculated MW	226.12
Formulation	C6H11O7P
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:O=C(O)C(CP(O)(O)=O)CCC(O)=O

Pathways:MetabolismReceptor:GCPIIBoling pt:Melting pt:Solubility:Water: 28 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Rais R, et al. J Pharm Biomed Anal. 2014 Jan;88:162-9.

Note: This product is for in vitro research use only