

## Vasopressin

Catalog No: #S3443


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Vasopressin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	1050.22
Formulation	C43H67N15O12S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular

systemSMILES:C[C@H](C)[C@H]1C(=O)N[C@H](C(=O)N[C@H](C(=O)N[C@@H](CSSC[C@@H](C(=O)N[C@H](C(=O)N1)CC2=CC=C(C(=C2)O)N)C(=O)N3CCC[C@H]3C(=O)N[C@@H](CCCN=C(N)N)C(=O)NCC(=O)N)CC(=O)N)CCC(=O)NPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in WaterAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Tsuji T, et al. J Physiol. 2017 Apr 12.

Note: This product is for in vitro research use only