

# Importazole


Catalog No: #S3445

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Importazole
Brief Description	Inhibitors
Purification	99.88%
Target Name	Others inhibitor
Calculated MW	318.42
Formulation	C <sub>20</sub> H <sub>22</sub> N <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CC(C)NC1=NC2=CC=CC=C2N=C1N3CCCC3  
 Pathways: Others  
 Receptor: importin-β  
 Boiling pt: Melting pt: Solubility: DMSO: 33 mg/ml  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Soderholm JF, et al. ACS Chem Biol. 2011 Jul 15;6(7):700-8.

Note: This product is for in vitro research use only