

Thiomyristoyl

Catalog No: #S3447


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thiomyristoyl
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sirtuin inhibitor
Calculated MW	581.85
Formulation	C34H51N3O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: O=C(OCC1=CC=CC=C1)N[C@@H](CCCCNC(CCCCCCCCCCCC)=S)C(NC2=CC=CC=C2)=O
 Pathways: Chromatin/Epigenetic;
 DNA Damage/DNA Repair
 Receptor: SIRT1; SIRT2
 Boiling pt: Melting pt: Solubility: Chloroform, DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jing H, et al. Cancer Cell. 2016 Mar 14;29(3):297-310.

Note: This product is for in vitro research use only