

## Bay 41-4109 (racemate)

Catalog No: #S3449


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bay 41-4109 (racemate)
Brief Description	Inhibitors
Purification	99.00%
Target Name	HBV inhibitor
Calculated MW	395.76
Formulation	C18H13ClF3N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: O=C(C1=C(C)NC(C2=NC=C(F)C=C2F)=NC1C3=CC=C(F)C=C3Cl)OC Pathways: Microbiology & virology Receptor: HBV Bolin  
 g pt: Melting pt: Solubility: DMSO: 28 mg/ml Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Weber O, et al. Antiviral Res. 2002 May;54(2):69-78.

Note: This product is for in vitro research use only