

HhAntag

Catalog No: #S3460


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HhAntag
Brief Description	Inhibitors
Purification	97.00%
Target Name	Hedgehog/Smoothened inhibitor
Calculated MW	450.923
Formulation	C24H23ClN4O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CN(C)C1=CC2=C(C=C1)N=C(N2)C3=C(C=CC(=C3)NC(=O)C4=CC(=CC(=C4)OC)OC)ClPathways:GPCR/G Protein;
Stem CellsReceptor:GliBoling pt:Melting pt:Solubility:DMSO10 mM Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Actis, Marcelo; Connelly, Michele C.; Mayasundari, Anand et al. A structure-activity relationship study of small-molecule inhibitors of GLI1-mediated transcription. Biopolymers (2011), 95(1), 24-30.

Note: This product is for in vitro research use only