

Vesnarinone

Catalog No: #S3465


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vesnarinone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	395.459
Formulation	C22H25N3O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:COC1=C(C=C(C=C1)C(=O)N2CCN(CC2)C3=CC4=C(C=C3)NC(=O)CC4)OC

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO10 mM Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Katayama Y, et al. Inhibitory effects of Vesnarinone on cloned cardiac delayed rectifier K(+) channels expressed in a mammalian cell line. J Pharmacol Exp Ther. 2000 Jul;294(1):339-46.

Note: This product is for in vitro research use only