

iPLA2 γ Polyclonal Antibody

Catalog No: #41072

Package Size: #41072-1 50ul #41072-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	iPLA2 γ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	iPLA2 γ Polyclonal Antibody detects endogenous levels of iPLA2 γ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human iPLA2 γ .
Target Name	iPLA2 γ
Other Names	PNPLA8; IPLA22; IPLA2G; BM-043; Calcium-independent phospholipase A2-gamma; Intracellular membrane-associated calcium-independent phospholipase A2 gamma; iPLA2-gamma; PNPLA-gamma; Patatin-like phospholipase domain-containing protein 8; iPLA
Accession No.	Swiss-Prot: Q9NP80NCBI Gene ID: 50640
Uniprot	Q9NP80
GeneID	50640;
SDS-PAGE MW	88kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images

Western Blot analysis of HeLa cells using iPLA2 γ Polyclonal Antibody

Note: This product is for in vitro research use only