


Description

Product Name	MD 39-AM
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	268.338
Formulation	C14H12N4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CSC1=NC2=C(C=CC=N2)C(=N1)NC3=CC=CC=C3

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: SolubleAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Prez MJ, et al. Preliminary pharmacokinetics of a new pyridopyrimidine derivative. *Arzneimittelforschung*. 1991 Jun;41(6):640-4.

Note: This product is for in vitro research use only