

Sudoterb free base

Catalog No: #S3480


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sudoterb free base
Brief Description	Inhibitors
Purification	98.00%
Target Name	ATPase inhibitor
Calculated MW	519.572
Formulation	C ₂₉ H ₂₈ F ₃ N ₅ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1=C(C=C(N1NC(=O)C2=CC=NC=C2)C3=CC=CC=C3)CN4CCN(CC4)C5=CC=CC(=C5)C(F)(F)F

Pathways: Membrane transporter/Ion channel Receptor: ATP Boiling pt: Melting pt: Solubility: CDCl₃ Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Didilescu C, Craiova UM. [Present and future in the use of anti-tubercular drugs]. *Pneumologia*. 2011 Oct-Dec;60(4):198-201.

Note: This product is for in vitro research use only