

## Tectorigenin

Catalog No: #S3488


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tectorigenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Apoptosis inducer inhibitor
Calculated MW	300.266
Formulation	C16H12O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:COC1=C(C=C2C(=C1O)C(=O)C(=CO2)C3=CC=C(C=C3)O)O

Pathways:ApoptosisReceptor:Apoptosis inducerBoiling pt:Melting pt:Solubility:Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Li QY, Chen L, Yan MM,et al.Tectorigenin regulates adipogenic differentiation and adipocytokines secretion via PPARγ and IKK/NF-κB signaling.Pharm Biol. 2015;53(11):1567-75.

Note: This product is for in vitro research use only