

Dasabuvir (ABT-333)

Catalog No: #S3489


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dasabuvir (ABT-333)
Brief Description	Inhibitors
Purification	98.00%
Target Name	HCV Protease inhibitor
Calculated MW	493.578
Formulation	C26H27N3O5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC(C)(C)C1=C(C(=CC(=C1)N2C=CC(=O)NC2=O)C3=CC4=C(C=C3)C=C(C=C4)NS(=O)(=O)C)OC

Pathways: Microbiology & virology; Proteases/Proteasome Receptor: NS5B polymerase

Bolting pt: Melting pt: Solubility: DMSO: 46 mg/mL
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Poordad F, et al. Exploratory study of oral combination antiviral therapy for hepatitis C. N Engl J Med. 2013 Jan 3;368(1):45-53.

Note: This product is for in vitro research use only