


## Description

Product Name	Danshensu
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	198.174
Formulation	C9H10O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:CC(C1=CC(=C(C=C1)O)O)(C(=O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Kwon G1, Kim HJ2, Park SJ1,et al.Anxiolytic-like effect of danshensu [(3-(3,4-dihydroxyphenyl)-lactic acid)] in mice.Life Sci. 2014 Apr 17;101(1-2):73-8.

Note: This product is for in vitro research use only