

Clofentezine

Catalog No: #S3493


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Clofentezine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	303.146
Formulation	C ₁₄ H ₈ Cl ₂ N ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C1=CC=C(C(=C1)C2=NN=C(N=N2)C3=CC=CC=C3Cl)Cl

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: < 8.2 mg/mL

Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Marcic D. The effects of clofentezine on life-table parameters in two-spotted spider mite Tetranychus urticae. Exp Appl Acarol. 2003;30(4):249-63.

Note: This product is for in vitro research use only