

Tasimelteon

Catalog No: #S3495


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tasimelteon
Brief Description	Inhibitors
Purification	99.00%
Target Name	Melatonin Receptor agonist
Calculated MW	245.322
Formulation	C15H19NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CCC(=O)NC[C@@H]1C[C@H]1C2=C3CCOC3=CC=C2
Pathways: GPCR/G Protein
Receptor: Melatonin Receptor
Boiling pt: Melting pt: Solubility: DMSO: 33 mg/mL

Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Leger D, et al. Safety profile of tasimelteon, a melatonin MT1 and MT2 receptor agonist: pooled safety analyses from six clinical studies. *Expert Opin Drug Saf.* 2015;14(11):1673-85.

Note: This product is for in vitro research use only