

## Remodelin

Catalog No: #S3499


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Remodelin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	282.365
Formulation	C15H14N4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:C1CCC(=NNC2=NC(=CS2)C3=CC=C(C=C3)C#N)C1

Pathways:OthersReceptor:NAT10 inhibitorBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Larrieu D, et al. Chemical inhibition of NAT10 corrects defects of laminopathic cells. Science. 2014 May 2;344(6183):527-32.

Note: This product is for in vitro research use only