

## MC-Val-Cit-PABC-PNP

Catalog No: #S3503


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	MC-Val-Cit-PABC-PNP
Brief Description	Inhibitors
Purification	95.00%
Target Name	AChR inhibitor
Calculated MW	737.767
Formulation	C35H43N7O11
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CC(C)[C@@H](C(=O)N[C@@H](CCCNC(=O)N)C(=O)NC1=CC=C(C=C1)COC(=O)OC2=CC=C(C=C2)[N+](=O)[O-]NC(=O)CC(C)CN3C(=O)C=CC3=O Pathways: Neuroscience Receptor: ADCs Linker Boling pt: Melting pt: Solubility: DMSO: 37 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Anti-cd70 antibody-drug conjugates and their use for the treatment of cancer and immune disorders

Note: This product is for in vitro research use only