

IRF-3 Polyclonal Antibody

Catalog No: #41075

Package Size: #41075-1 50ul #41075-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRF-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	IRF-3 Polyclonal Antibody detects endogenous levels of IRF-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human IRF-3 around the non-phosphorylation site of S396.
Target Name	IRF-3
Other Names	IRF3; Interferon regulatory factor 3; IRF-3
Accession No.	Swiss-Prot: Q14653NCBI Gene ID: 3661
Uniprot	Q14653
GeneID	3661;
SDS-PAGE MW	47kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using IRF-3 Polyclonal Antibody

Note: This product is for in vitro research use only