

tert-Butyl ((3-hydroxyazetidin-3-yl)methyl)carbamate

Catalog No: #S3505


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	tert-Butyl ((3-hydroxyazetidin-3-yl)methyl)carbamate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	202.254
Formulation	C ₉ H ₁₈ N ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC(C)(C)OC(=O)NCC1(CNC1)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:

SolubleAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only