

2-CHLORO-5-CYANOPYRAZINE

Catalog No: #S3507


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-CHLORO-5-CYANOPYRAZINE
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	139.542
Formulation	C5H2ClN3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=C(N=CC(=N1)Cl)C#N Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO:

Soluble Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Laixing Hu, et al. Synthesis and activity of azaterphenyl diamidines against *Trypanosoma brucei rhodesiense* and *Plasmodium falciparum*. *Bioorg Med Chem*. Author manuscript; available in PMC 2013 Oct 30.

Note: This product is for in vitro research use only